

INOGEN *at* **HOME**[®]
Home Oxygen Concentrator



Inogen At Home System (GS-100):

Only 18 Pounds
5 Liters Per Minute, Continuous Flow
Energy Efficient

INOGEN *at* HOME[®]

The Inogen At Home is one of the lightest and energy-efficient 5 liter per minute continuous flow oxygen concentrators on the market today.

At approximately 18 pounds, it is half the weight of some other oxygen concentrators in use today. At low power consumption levels, the Inogen At Home is an energy-efficient choice for your next home oxygen concentrator.

Weight	18 lbs.
Size	Width: 13.0" Depth: 7.0" Height: 16.5"
Oxygen Flow	Continuous flow system Five flow settings from 1-5 liters per minute
Power Sensing	AC Power Supply 100-240V, 275W, 50-60Hz (auto sensing to allow worldwide use)
Warranty	3-year warranty*
Operation	Simple control functions
Use	Designed for stationary use

*See Inogen At Home Limited Warranty Statement for complete warranty information

Manufactured by:

inogen[®]

301 Coromar address
Goleta, California 93117
www.inogen.com

