



OXYGEN CONCENTRATORS

Doctor Discussion Guide



Call 1-800-819-0646
for a free consultation

INOGEN®
PORTABLE OXYGEN CONCENTRATORS

Oxygen therapy that moves with you.

Introducing the Inogen Rove Family

If you require and have been prescribed supplemental oxygen therapy, get ready for a welcome breath of fresh air with the Inogen Rove 4™ and Rove 6™.

For over 20 years, Inogen has helped over 1 million patients replace bulky portable oxygen tanks with lightweight, portable oxygen concentrators (POCs) that deliver supplemental oxygen while on the go.

As America's leader in portable oxygen therapy, Inogen helps patients receive the benefits of oxygen therapy while maintaining their active lifestyles. This guide will introduce you to our line of oxygen therapy solutions. **There's only *one* Inogen.**

Call 1-800-819-0646

to speak with an oxygen specialist
and confirm the right option for you.



INOGEN®
PORTABLE OXYGEN CONCENTRATORS

Our compact alternative to tanks.

Portable oxygen concentrators (POCs) provide a safe source of oxygen-enriched air to deliver medical-grade oxygen as prescribed by a physician, freeing patients from managing heavy, cumbersome tanks and planning for refills.

Inogen POCs work by pulling in the air around you, separating the nitrogen, and delivering a concentrated amount of oxygen on each breath.

With reliable, pulse-dose based Intelligent Delivery Technology®, Inogen POCs are powered to provide medical grade oxygen delivery all day, every day. Pulse dose oxygen is delivered based primarily on breathing rate as a sensor determines when a patient's inhale begins to deliver an oxygen pulse at that moment. Patients get the oxygen they need, when they need it.

- 1 At flow setting 1
- 2 See limited warranty statement on www.inogen.com for complete warranty information.
- 3 At flow setting 2
- 4 When used with a standard battery

Call **1-800-819-0646** to get started

Inogen Rove 4™



7.50 in.

- Only 2.9 lbs. (with standard battery)
- 4 Pulse-dose, Intelligent Delivery Technology® Settings (1-4)
- Up to 3.00 hours run time with standard battery; up to 5.75 hours with extended battery¹
- 3-year or lifetime warranty²
- 39 dBA³ operating volume
- AC Power Supply 100-240 V, 100 W, 50-60 Hz; Auto sensing to allow for worldwide use
- 5.9" L x 2.7" W

Inogen Rove 6™



8.14 in.

- Only 4.8 lbs (with standard battery)
- 6 Pulse-dose, Intelligent Delivery Technology® Settings (1-6)
- Up to 6.25 hours run time with standard battery; up to 12.75 hours with extended battery¹
- 3-year or lifetime warranty²
- 39 dBA³ operating volume
- AC Power Supply 100-240 V, 100 W, 50-60 Hz; Auto sensing to allow for worldwide use
- 7.18" L x 3.27" W

INOGEN®
PORTABLE OXYGEN CONCENTRATORS

Accessories for life on-the-go.

Purchase carry bags, extra batteries, power cables, replacement filters, and more accessories to keep you moving at inogen.com.

Rove 6 Carry Bag



Rove 4 Extended Batteries



All Inogen portable systems meet FAA acceptance criteria for carriage and use onboard aircraft.*

Call **1-800-819-0646** to get started

* Based on technical performance test data.


“ My Inogen system is so much lighter than canisters, and gives me hours of freedom to do what I want.

Jerry L.**

Go hands free
with the backpack



** The opinions expressed herein are independently prepared and for informational purposes only. The views expressed are of the individuals' only. The individual featured in this content did not receive any financial compensation for providing their testimonial/s to Inogen, Inc. to use in our program.

A person wearing an orange long-sleeved shirt is walking away from the camera on a path in a sunlit park. The background is filled with out-of-focus green trees and foliage, creating a bokeh effect. The lighting is bright and warm, suggesting a sunny day. On the left side of the image, there are several overlapping translucent blue and white geometric shapes that form a decorative border.

“ I’m beyond satisfied with my Inogen[®] portable and plug-in concentrators. The portable one lets me keep doing the work I enjoy, and walks and light exercise, while the plug-in is so quiet and light weight. It travels with me to visit family, and on trips with my husband. I’m so thankful.

Debra E.*

* The opinions expressed herein are independently prepared and for informational purposes only. The views expressed are of the individuals' only. The individual featured in this content did not receive any financial compensation for providing their testimonial/s to Inogen, Inc. to use in our program.





INOGEN®
STATIONARY OXYGEN CONCENTRATOR

Compact and quiet at-home therapy.

Inogen At Home®

Our Inogen At Home system combines 5-liter per minute continuous flow oxygen delivery with technology that is quiet, letting you watch TV or visit with friends with little noise disturbance.


At just 18 pounds, Inogen At Home is a compact unit you can keep at your bedside, pack in carry on luggage, or use it during the day as a companion around the house. The At Home is energy-efficient to help you save in electric costs over time.

Call 1-800-819-0646

to speak with an oxygen specialist
and confirm the right option for you.

Inogen At Home[®]

- 18.0 lbs.
- 5 Continuous flow settings (1-5)
- 3-year or lifetime warranty¹
- 45 dBA² operating volume
- AC Power Supply 100-240 V, 275 W, 50-60 Hz;
Auto sensing to allow for worldwide use



INOGEN®
STATIONARY OXYGEN CONCENTRATOR

Power & efficiency in a compact size.

Save on your electric bill each year*

Thanks to its energy technology, Inogen At Home® is an energy efficient option for oxygen concentrators used at home. Because it's designed for 24/7 use, the power consumption can significantly add up.



* Calculations are based on published US average hourly kilowatt rates, and assume stationary oxygen concentrator usage of 20 hours per day. A flow setting of 2 liters per minute was used as the Inogen At Home® power draw, with the published specifications of a common competitive stationary concentrator used as the standard market concentrator power draw. Actual costs and/or savings may vary.

1 See limited warranty statement on www.inogen.com for complete warranty information.

2 At flow setting 2

Call **1-800-819-0646** to get started



Joining the Inogen Family



Inogen's 5-Star Affordable Access Options

We're committed to providing affordable oxygen therapy for patients, that's why our oxygen specialists will guide you through our 5-Star Affordable Access Options.

- **Medicare:** A POC may be covered by Medicare as an alternative to tanks for portable oxygen therapy.
- **Private Insurance:** A POC may be covered by your private insurance for your portable oxygen therapy needs.
- **HSA/FSA Funds:** If an Inogen customer has been determined out of network, HSA/FSA funds may be used to purchase directly.
- **Financing:** You may be eligible for financing to purchase an Inogen POC of your choice.
- **Cash Purchase:** Keep your provider and receive fast shipping for a POC of your choice shipped directly to you.



Next Steps

- Your doctor provides a prescription saying that you require supplemental oxygen and/or have a severe lung disease.
- Your medical documentation indicates that you are mobile in your home and would benefit from the use of a portable system.
- Skip the tank. Inogen concentrators are available day one of your oxygen therapy when prescribed a POC.
- Call **1-800-940-6125**. We'll check your portable oxygen concentrator eligibility with Medicare, private insurance, or walk you through cash pay options.



SCAN THE CODE

to discover how people like
you are enjoying their mobility
with Inogen®



Inogen, Inc.
859 Ward Drive, Suite 200, Goleta, California 93111
Copyright © 2025 Inogen | EX_EN_USA | PM250036

1-800-819-0646
Call for a free consultation